

## Guide for Reading

### Key Concepts

- What are the products of the reaction of an acid with a base?
- What is the endpoint of a titration?

### Vocabulary

neutralization reactions  
equivalence point  
titration  
standard solution  
end point

### Reading Strategy

#### Identifying a Sequence

A sequence is the order in which a series of events occurs. As you read about acid-base titrations, list the steps that should be used in carrying out a precise titration. Include the reactants and how they are measured, the indicator, and what to look for as the titration nears its end point.

**Figure 19.19** For fish and other aquatic animals to survive, the water in which they live must be maintained at the proper pH.


## Connecting to Your World

Nearly all of the adult population suffers from acid indigestion at some time. Although hydrochloric acid is always present in the stomach, an excess can cause heartburn and a feeling of nausea. A common way to relieve the pain of acid indigestion is to take antacids to neutralize the stomach acid. The active ingredient in many antacids is sodium hydrogen carbonate, aluminum hydroxide, or magnesium hydroxide. In this section, you will learn what a neutralization reaction is.



## Acid-Base Reactions

If you mix a solution of a strong acid containing hydronium (hydrogen ions) with a solution of a strong base that has an equal number of hydroxide ions, a neutral solution results. The final solution has properties that are characteristic of neither an acidic nor a basic solution. Consider these examples:



In each example, a strong acid reacts with a strong base. If solutions of these substances are mixed in the mole ratios specified by the balanced equation, neutral solutions will result. Similar reactions of weak acids and/or weak bases do not usually produce neutral solutions. In general, however, reactions in which an acid and a base react in an aqueous solution to produce a salt and water are called **neutralization reactions**. The formation of water in a neutralization reaction is shown in Figure 19.20.


Neutralization reactions are one way to prepare pure samples of salts. You could prepare potassium chloride, for example, by mixing equal molar quantities of hydrochloric acid and potassium hydroxide. An aqueous solution of potassium chloride would result. You could heat the solution to evaporate the water, leaving the salt potassium chloride. Table 19.9 lists some common salts and their applications.

**In general, the reaction of an acid with a base produces water and one of a class of compounds called salts.** When you hear the word *salt*, you may think of the substance that flavors your French fries or scrambled eggs.

Table salt (sodium chloride) is one example of a salt, but there are many more. Salts are compounds consisting of an anion from an acid and a cation from a base.



Name	Formula	Applications
Ammonium sulfate	$(\text{NH}_4)_2\text{SO}_4$	Fertilizer
Barium sulfate	$\text{BaSO}_4$	Gastrointestinal studies; white pigment
Calcium chloride	$\text{CaCl}_2$	De-icing roadways and sidewalks
Potassium chloride	$\text{KCl}$	Sodium-free salt substitute
Silver bromide	$\text{AgBr}$	Photographic emulsions
Sodium hydrogen carbonate (baking soda)	$\text{NaHCO}_3$	Antacid
Sodium carbonate decahydrate (washing soda)	$\text{Na}_2\text{CO}_3 \cdot 10\text{H}_2\text{O}$	Glass manufacture; water softener
Sodium chloride (table salt)	$\text{NaCl}$	Body electrolyte; chlorine manufacture




**Figure 19.20** In a neutralization reaction, hydronium ions ( $\text{H}_3\text{O}^+$ ) combine with hydroxide ions ( $\text{OH}^-$ ) to form neutral water.


The properties of acids, bases, and salts help explain many diverse phenomena. The usefulness of antacids, for example, depends on the process of acid-base neutralization. Farmers use a similar process to control the pH of soil.

## Titration


Acids and bases sometimes, but not always, react in a 1:1 mole ratio.



When sulfuric acid reacts with sodium hydroxide, however, the ratio is 1:2. Two moles of the base sodium hydroxide are required to neutralize one mole of  $\text{H}_2\text{SO}_4$ .



Similarly, hydrochloric acid and calcium hydroxide react in a 2:1 ratio.



You will notice in the preceding examples that the number of moles of hydrogen ions provided by the acid are equivalent to the number of hydroxide ions provided by the base. When an acid and base are mixed, the **equivalence point** is when the number of moles of hydrogen ions equals the number of moles of hydroxide ions.



### Checkpoint

**What is the equivalence point of a reaction between an acid and a base?**

How many moles of sulfuric acid are required to neutralize 0.50 mol of sodium hydroxide?

**1 Analyze** *List the knowns and the unknown.*

**Knowns**

- mol NaOH = 0.50 mol
- $\text{H}_2\text{SO}_4(\text{aq}) + 2\text{NaOH}(\text{aq}) \longrightarrow \text{Na}_2\text{SO}_4(\text{aq}) + 2\text{H}_2\text{O}(\text{l})$
- $\frac{\text{mol H}_2\text{SO}_4}{\text{mol NaOH}} = \frac{1}{2}$

**Unknown**

- moles  $\text{H}_2\text{SO}_4 = ?$  mol

**2 Calculate** *Solve for the unknown.*

The mole ratio of  $\text{H}_2\text{SO}_4$  to NaOH is 1:2. The necessary number of moles of  $\text{H}_2\text{SO}_4$  is calculated using this ratio.

$$0.50 \text{ mol NaOH} \times \frac{1 \text{ mol H}_2\text{SO}_4}{2 \text{ mol NaOH}} = 0.25 \text{ mol H}_2\text{SO}_4$$

**3 Evaluate** *Does the result make sense?*

Because the mole ratio of  $\text{H}_2\text{SO}_4$  to NaOH is 1:2, the expected number of moles of  $\text{H}_2\text{SO}_4$  should be half the given number of moles of NaOH. The answer should have two significant figures.

**Practice Problems**

**30.** How many moles of potassium hydroxide are needed to completely neutralize 1.56 mol of phosphoric acid?

**31.** How many moles of sodium hydroxide are required to neutralize 0.20 mol of nitric acid?

**Math**

**Handbook**

For help with dimensional analysis, go to page R66.

**Interactive  
Textbook**


**Problem-Solving 19.30** Solve Problem 30 with the help of an interactive guided tutorial.

with **ChemASAP**

You can determine the concentration of acid (or base) in a solution by performing a neutralization reaction. You must use an appropriate acid-base indicator to show when neutralization has occurred. As you can see in Figure 19.21, the juice of the red cabbage is an acid-base indicator. In the laboratory, phenolphthalein is often the preferred indicator for acid-base neutralization reactions. Solutions that contain phenolphthalein turn from colorless to deep pink as the pH of the solution changes from acidic to basic. In slightly basic solutions, the indicator is very faintly pink.

**Figure 19.21** Red cabbage juice is used as an acid-base indicator. As the solution changes from highly acidic to basic, the color changes from red to violet to green to yellow. **Predicting** *Would the yellow solution have a high or low pH?*






acid (plus a few drops of phenolphthalein indicator) in a flask is placed beneath a buret filled with a base of known concentration. **(b)** Base is slowly added from the buret to the acid while the flask is gently swirled. **(c)** A change in the color of the indicator signals that neutralization has occurred.

The steps in a neutralization reaction are as follows.

1. A measured volume of an acid solution of unknown concentration is added to a flask.
2. Several drops of the indicator are added to the solution while the flask is gently swirled.
3. Measured volumes of a base of known concentration are mixed into the acid until the indicator just barely changes color.

The process of adding a known amount of solution of known concentration to determine the concentration of another solution is called **titration**. The solution of known concentration is called the **standard solution**. Use a buret to add the standard solution. A titration is continued until the indicator shows that neutralization has just occurred. The point at which the indicator changes color is the **end point** of the titration. The titration of an acid of unknown concentration with a standard base is shown in Figure 19.22. You can use a similar procedure to find the concentration of a base using a standard acid.

Figure 19.23 shows how the pH of a solution changes during the titration of a strong acid (HCl) with a strong base (NaOH). The pH of the initial acid solution is low. As the base is added, the pH increases because some of the acid is neutralized. As the titration approaches the point of neutralization, at a pH of 7, the pH increases dramatically as hydrogen ions are used up. Once past the point of neutralization, additional base produces a further increase of pH. **Key** The point of neutralization is the **end point of the titration**. At this point, the contents of the beaker consist of only  $\text{H}_2\text{O}$  and NaCl, which is the salt derived from the strong acid HCl and the strong base NaOH, plus a trace of indicator.



**Figure 19.23** In this titration of a strong acid with a strong base, 0.10M NaOH is slowly added from a buret to 50.0 mL of 0.10M HCl in the beaker. The equivalence point, the midpoint on the vertical portion of the pH titration curve, occurs at 50.0 mL of NaOH added. **Interpreting Diagrams** What is true concerning  $[\text{H}^+]$  and  $[\text{OH}^-]$  at the equivalence point?

### Interactive Textbook

**Simulation 26** Simulate the titration of several acids and bases and observe patterns in the pH at equivalence.

with **ChemASAP**



A 25-mL solution of  $\text{H}_2\text{SO}_4$  is completely neutralized by 18 mL of 1.0M NaOH. What is the concentration of the  $\text{H}_2\text{SO}_4$  solution?

**1 Analyze** *List the knowns and the unknown.*

**Knowns**

- molarity base = 1.0M NaOH
- volume base = 18 mL = 0.018 L
- volume acid = 25 mL = 0.025 L
- $\text{H}_2\text{SO}_4(aq) + 2\text{NaOH}(aq) \longrightarrow \text{Na}_2\text{SO}_4(aq) + 2\text{H}_2\text{O}(l)$

**Unknown**

- molarity acid = ?M  $\text{H}_2\text{SO}_4$

Use the molarity to convert the volume of base to moles of base. Use the mole ratio to find the number of moles of  $\text{H}_2\text{SO}_4$ . Calculate the molarity by dividing the number of moles of  $\text{H}_2\text{SO}_4$  by the number of liters of  $\text{H}_2\text{SO}_4$ .

The conversion steps are as follows:



**2 Calculate** *Solve for the unknown.*

$$0.018 \text{ L NaOH} \times \frac{1.0 \text{ mol NaOH}}{1 \text{ L NaOH}} \times \frac{1 \text{ mol H}_2\text{SO}_4}{2 \text{ mol NaOH}} = 0.0090 \text{ mol H}_2\text{SO}_4$$

$$\text{molarity} = \frac{\text{moles}}{\text{liters}} = \frac{0.0090 \text{ mol}}{0.025 \text{ L}} = 0.36 \text{ M}$$

The  $[\text{H}_2\text{SO}_4]$  is 0.36M.

**3 Evaluate** *Does the result make sense?*

Because the volume of acid was greater than the volume of base, the concentration is less than 1.0M. The answer has two significant figures.

**Practice Problems**

32. How many milliliters of 0.45M HCl will neutralize 25.0 mL of 1.00M KOH?

33. What is the molarity of  $\text{H}_3\text{PO}_4$  if 15.0 mL is completely neutralized by 38.5 mL of 0.150M NaOH?

**Math**

**Handbook**

For help with dimensional analysis, go to page R66.



**Problem-Solving 19.33**

Solve Problem 33 with the help of an interactive guided tutorial.

with **ChemASAP**

## 19.4 Section Assessment

- Key Concept** What are the products of a reaction between an acid and a base?
- Key Concept** What occurs at the endpoint of a titration?
- How many moles of HCl are required to neutralize aqueous solutions of these bases?
  - 2 mol  $\text{NH}_3$
  - 0.1 mol  $\text{Ca}(\text{OH})_2$
- Write complete balanced equations for the following acid-base reactions?
  - $\text{H}_2\text{SO}_4(aq) + \text{KOH}(aq) \longrightarrow$
  - $\text{H}_3\text{PO}_4(aq) + \text{Ca}(\text{OH})_2(aq) \longrightarrow$
  - $\text{HNO}_3(aq) + \text{Mg}(\text{OH})_2(aq) \longrightarrow$

**Connecting**

**Concepts**

**Types of Reactions** Reread the information on types of chemical reactions in Section 11.2. Which type of reaction are all acid-base neutralizations? Explain your answer.



**Assessment 19.4** Test yourself on the concepts in Section 19.4.

with **ChemASAP**